



MOTOR: THREE SPEED FORWARD CURVE MOTOR/BLOWER
 1/5 HP 115V, 60 HZ, 1 PHASE 3.0 F.L.A
 1/5 HP 230V, 50/60 HZ, 1 PHASE 1.5 F.L.A
 1/5 HP 277V, 60 HZ, 1 PHASE 1.2 F.L.A

FILTER: HEPA FILTER MEDIA 99.99% @ 0.3 MICRON
 ULPA FILTER MEDIA 99.9995% @ 0.12 MICRON
 UL LISTED: UL 507, 2F82 AIR FILTERING APPLIANCES, E152685



5601 BALLOON FIESTA PKWY N.E.
 ALBUQUERQUE NM 87113
 PHONE: 1 . 800 . 545 . 6598
 FAX: 1 . 505 . 345 . 8875
 www.envirco.com
 email: envirco@envirco.com

MAC10 ORIGINAL 2X2 STD

UNLESS SPECIFIED,
 ALL DIMENSIONS
 ARE IN INCHES.
 TOLERANCE
 X .XX .XXX
 1 .03 .010
 UNLESS NOTED

THIS DRAWING IS THE PROPERTY
 OF ENVIRCO CORP., AND MAY NOT
 BE REPRODUCED OR PUBLISHED
 WITHOUT PERMISSION FROM
 ENVIRCO CORP.

COPYRIGHT 2001 ENVIRCO CORP.

DRAWING NO: XL 2x2	
DATE: 11/28/01	DR: MWM
REV: A	BY: MWM
DESIGN PLAN:	REV DATE: 11/28/01
	N/A